

NELLA™
VULIGHT®

Single-Use LED Illuminator
Developed with OB-GYN Providers and Patients



The Small and Mighty Speculum Light



- Shadow-free illumination
- Fits most specula, metal or plastic
- Rotates along top bill – not blocked by tools
- Reduced risk of cross contamination compared to gooseneck light

44% Brighter

the Nella VuLight has been shown to have 44% higher Lux than a gooseneck exam light¹



Shadow-free:
Light where you want it.



Targeted Illumination
Nella VuLight
Single Use LED Light



Traditional Illumination
Welch Allyn 48810
Exam Light

- **Brighter illumination** compared to a standard gooseneck light¹
- **Close to cervix** placement provides direct illumination even in patients with difficult anatomy, excess tissue or first-time patients
- **Securely reposition** along top bill

Use VuLight for:

General gynecology exams
LEEP procedures
Biopsies
IUD insertions/removals
Other in-office procedures
Emergency department

Provider Needs

Provider Needs	VuLight	Gooseneck exam light
Focused cervical illumination ¹	✓	✗
Higher Lux (brightness) ¹	✓	✗
Shadow-free ¹	✓	✗
No risk of cross-contamination	✓	✗
No need to reposition external light source	✓	✗
Portable solution	✓	✗
No provider contortion necessary	✓	✗

¹ Ceek Women's Health (2020). VuLight Comparison Verification. Internal Report: unpublished.



Products for a better pelvic exam.



reddot winner 2020



Ordering Information:

VL050 Nella VuLight Single-Use LED Illuminator (Box of 50)

order@ceekwomenshealth.com or call **1-855-559-CEEK**

1455 NW Irving Street, Suite 200, Portland, Oregon 97209