

NELLA™
NUSPEC™

Reusable Vaginal Speculum
Developed with OB-GYN Providers and Patients



The Smarter Reusable Speculum



- Improved visibility with integrated sidewall retractors
- Temperature neutral- no need to warm
- Quiet operation keeps patients at ease
- Smooth- no lubricant required
- Strong, autoclavable polymer

99%

of women prefer the NuSpec
over a metal speculum¹

Full Visibility Meets Ultimate Patient Comfort



NuSpec



sidewall tissue retracted

Pederson



sidewall tissue encroachment

Graves



sidewall tissue encroachment

Patient: 41 y.o., BMI 30, g2p2

Patient and Provider Preferred.

- **Improved cervical visibility**² for patients with sidewall laxity
- **Increased comfort**² for all patients, especially adolescent, menopausal or cancer survivors
- **Minimized anxiety and discomfort**² for first-time or sensitive patients

Use NuSpec for:

General gynecology exams
Colposcopy
Biopsies
IUD insertion & removal
Other in-office procedures

Provider Needs

	NuSpec	Metal
Sidewall retention	✓	✗
Easy adjustment in use	✓	✗
Ergonomic handle	✓	✗
One-handed locking	✓	✗
No disassembly required	✓	✗

Patient Needs

	NuSpec	Metal
Smaller bill size	✓	✗
Quiet operation	✓	✗
Neutral temperature	✓	✗
Smooth material (no lubricant required)	✓	✗

¹ Ceek Women's Health (2020). *NuSpec Consumer Preference Study*. Internal Report.

² Ceek Women's Health (2017). *VIOLET Clinical Trial*. ClinicalTrials.gov Identifier: NCT03123367. Internal Report.



Products for a better pelvic exam.

Ordering Information:

NU010 Nella NuSpec Reusable Vaginal Speculum (Box of 10)

NU001 Nella NuSpec Reusable Vaginal Speculum (Single Speculum)

order@ceekwomenshealth.com or call **1-855-559-CEEK**

1455 NW Irving Street, Suite 200, Portland, Oregon 97209

© 2018-2020 Ceek Women's Health NUUSM1R3-2020