



**Your partner  
for medical furniture.**



## A Sense of Company

With more than 90 years of experience, SCHMITZ is one of the leading international manufacturers of medical products.

A family-owned business with a living tradition. Our typical employee has been with the company for some 18 years, proof of how strongly they identify with our values and philosophy.

In 2012, the company's management passed to Friedrich and Matthias Schmitz, the fourth generation of the family.

We are proud to be located in Germany and are deeply rooted in our home region. We develop and manufacture here.

90% of our suppliers are from Germany, of which 60% are from the local region.

## Made in Germany



---

*'The transition to the next generation is always something special.'*





## **A Sense of Care**

We understand the services which are our lifeblood, and aim to achieve maximum customer satisfaction. From consulting through planning to installation.


Cost-effectiveness through longevity is what you can expect from us.

With a well thought-out service concept and comprehensive support, SCHMITZ works alongside you as partner from start to finish.

### **Technical Customer Service**

T +49 (0)2377 84 549

service@schmitz-soehne.de



- Over 80 international trade partners
- Sales offices

## SCHMITZ – represented worldwide









**varimed®**

**Couches, stools, multipurpose trolleys  
and examination seats**

**The indispensable helper in your daily routine.**

**Quality and longevity**

High-quality materials combined with outstanding workmanship guarantee the longevity of our products.

**The everyday helper**

medi-net multipurpose trolley – our reliable helper that you can count on throughout your workday.







---

Examination couch 211.1080.0  
upholstery gold orange 106, powder coating grey aluminium RAL 9007



---

ECG examination couch 211.1310.0  
upholstery lavender purple 107, powder coating white aluminium RAL 9006





Examination and massage couch 211.9301.0  
upholstery carmine red 92, powder coating grey aluminium RAL 9007



Recovery couch 211.9501.0  
upholstery blossom white 84, frame chromium-plated




Examination couch 211.1180.0  
upholstery ice blue 108, powder coating grey white RAL 9002  
Examination couch 211.1190.0  
upholstery agave green 109, frame chromium-plated

## varimed® Couches


<sup>1</sup> classic  
<sup>2</sup> comfort

Model No. Description



### 211.1080.0 Examination couch<sup>123</sup>

rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, 2-section 50 mm thick foam rubber upholstery, levelling feet, head section adjustable up to +40°. Safe working load: 225 kg.



### 211.1180.0 Examination couch<sup>123</sup>

rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, 2-section 70 mm thick foam rubber upholstery, levelling feet, head section adjustable up to +40°. Safe working load: 225 kg.

### 211.1190.0 Examination couch<sup>12</sup>

chromium-plated frame, otherwise like 211.1180.0



### 200.0970.0 Paper roll holder

for paper roll up to 500 mm width



### 200.1030.0 Paper roll

500 mm wide



### 200.6500.0 Coloured frame

(w/o ill.) acc. to colour chart, for 211.1080.0 and 211.1180.0

**Surcharge**

#### Note:

For orthopaedic and massage treatment, couches with transverse and longitudinal bracings are recommended!

**All examination couches with spacer block between frame and upholstery!**

### 200.4380.0 Head section with nose inlet


(w/o ill.) for 211.1180.0 and 211.1190.0

**Surcharge**


### 200.4390.0 Filling for nose inlet

(w/o ill.)

- <sup>1</sup> 4 colours for leatherette upholstery in "classic" quality without surcharge.
- <sup>2</sup> 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.
- <sup>3</sup> Standard powder-coating colours for frame: grey-white 9002, white aluminium RAL 9006 or grey aluminium RAL 9007. Please indicate desired colour upon order placement.

varimed® **Couches**<sup>1</sup> classic  
<sup>2</sup> comfort

Model No. Description

**211.1310.0 ECG Examination couch<sup>123</sup>**   
800 mm wide, rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, with lateral bracing, 70 mm thick foam rubber upholstery with 2-section cover, levelling feet, head section adjustable up to +40°. Safe working load: 225 kg.

**211.1320.0 ECG Examination couch<sup>12</sup>**   
chromium-plated frame, otherwise like 211.1310.0



**200.0940.0 Paper roll holder**  
for 800 mm wide couches



**200.1030.0 Paper roll**  
500 mm wide

**200.5100.0 Longitudinal bracing**  
(w/o ill.) powder coated RAL 9002, RAL 9006 or RAL 9007,  
for 211.1310.0

**200.5110.0 Longitudinal bracing**  
(w/o ill.) chromium-plated, for 211.1320.0

**200.5120.0 Longitudinal bracing**  
(w/o ill.) powder coated according to colour chart

**200.6500.0 Coloured frame**  
(w/o ill.) acc. to colour chart, for 211.1310.0

**Surcharge**

<sup>1</sup> 4 colours for leatherette upholstery in “classic” quality without surcharge.

<sup>2</sup> 15 colours for leatherette upholstery in “comfort” quality (surcharge). Please see page 40.

<sup>3</sup> Standard powder-coating colours for frame: grey-white 9002, white aluminium RAL 9006 or grey aluminium RAL 9007.  
Please indicate desired colour upon order placement.



Model No.	Description	
<b>200.0000.0</b> (w/o ill.)	<b>Different height</b> up to max. 830 mm (from 700 mm height upwards, including transverse and longitudinal bracing) for couches 211.1180.0 – 211.1320.0	<b>Surcharge</b>
<b>200.0010.0</b> (w/o ill.)	<b>Mobility</b> for couches 211.1180.0 - 211.1320.0, 4 swivel-mounted castors diam. 125 mm, 2 castors with total locking device, antistatic, 1 castor with directional locking, with revolving bumpers diam. 150 mm, longitudinal bracing. All castors are resistant to abrasions. (please indicate desired final height of couch!)	<b>Surcharge</b>
<b>200.0020.0</b> (w/o ill.)	<b>Different width<sup>4</sup></b> of +/- 150 mm (max. width of 900 mm) for couches 211.1180.0 – 211.1320.0	<b>Surcharge<sup>4</sup></b>
<b>200.0030.0</b> (w/o ill.)	<b>Different length<sup>4</sup></b> of +/- 150 mm for couches 211.1180.0 – 211.1320.0	<b>Surcharge<sup>4</sup></b>

<sup>4</sup> If both special height and width are required,  
the surcharge will only be applicable once.




Bariatric couch 211.1330.0  
upholstery caribbean blue 82, frame chromium-plated

## varimed® Couches

<sup>1</sup> classic  
<sup>2</sup> comfort

Model No. Description



### 211.1330.0 Bariatric couch<sup>12</sup>

800 mm wide, chromium-plated rectangular steel tube frame, with transverse and longitudinal bracing, 2-section mattress, 70 mm thick foam rubber upholstery, levelling feet, head section adjustable up to +40°  
Safe working load: 300 kg. Mobile version and special widths/lengths not available!



200.0940.0 Paper roll holder  
for 800 mm wide couches




200.1030.0 Paper roll  
500 mm wide

- <sup>1</sup> 4 colours for leatherette upholstery in “classic” quality without surcharge.  
<sup>2</sup> 15 colours for leatherette upholstery in “comfort” quality (surcharge). Please see page 40.

varimed® CEMD

**Mobile couches, recovery couches**<sup>1</sup> classic  
<sup>2</sup> comfort


Model No. Description

**211.9301.0 Examination and massage couch<sup>123</sup>**

with transverse and longitudinal bracing, rectangular steel tube frame, powder coated RAL 9002, RAL 9006 or RAL 9007, with 4 swivel-mounted castors diam. 125 mm, 2 castors with total locking device, antistatic, 1 castor with directional locking, with revolving bumpers diam. 150 mm, 70 mm foam rubber upholstery with 2-section cover, head section adjustable up to +40°. Safe working load: 225 kg.

**211.9351.0 Examination and massage couch<sup>12</sup>**

chromium-plated frame, otherwise like 211.9301.0.

**211.9501.0 Recovery couch<sup>12</sup>**

chromium-plated rectangular steel tube frame, with transverse and longitudinal bracing, with 4 swivel castors diam. 125 mm, 2 castors with total locking device, antistatic, 1 castor with directional locking, with revolving bumpers diam. 150 mm. 70 mm thick foam rubber upholstery with 2-section cover, head section adjustable up to +40°. Equipped with 2 side guards with attachments, chromium-plated, can be lowered down to height of upholstery. Safe working load: 225 kg.

**211.9511.0 Recovery couch<sup>12</sup>**

chromium-plated frame, 680 mm height, otherwise like 211.9501.0.

**All examination couches with spacer block between frame and upholstery!**

<sup>1</sup> 4 colours for leatherette upholstery in "classic" quality without surcharge.

<sup>2</sup> 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.

<sup>3</sup> Standard powder-coating colours for frame: grey-white 9002, white aluminium RAL 9006 or grey aluminium RAL 9007. Please indicate desired colour upon order placement.






Model No. Description


**200.0970.0 Paper roll holder**  
for paper roll up to 500 mm width.



**200.1030.0 Paper roll**  
500 mm wide.




**200.4900.0 1 Side guard and 1 side rail piece**  
side guard laterally adjustable on side rail,  
chromium-plated, with clamp, fitted to  
rectangular steel tube frame couches,  
may be lowered down to height of upholstery. **piece**



**200.4950.0 1 Side guard with clamp**  
screwed, rigid connection, fitted to rectangular  
steel tube frame couches, chromium-plated,  
may be lowered down to height of upholstery. **piece**

**200.0020.0 Different width of 800 mm**  
(w/o ill.) for couches 211.9301.0 to 211.9511.0. **Surcharge**



**200.5500.0 Storage basket**  
formed to nest onto transverse or longitudinal  
bracing, basket 520 × 490 × 140 mm,  
chromium-plated.

**200.6500.0 Coloured frame**  
(w/o ill.) acc. to colour chart, for 211.9301.0. **Surcharge**

All castors are non-marking and resistant to abrasions.



- Large and ergonomic height adjustment range
- No displacement of table-top during height adjustment
- Moulded base provides generous legspace




varimed® C€MD

## Height adjustable examination and treatment couches

<sup>1</sup> classic  
<sup>2</sup> comfort

Model No. Description




### 211.4100.0 Examination and treatment couch<sup>123</sup>

hydraulic height adjustment 500 to 920 mm by means of foot lever. Servo-assisted head section adjustment from  $-20^{\circ}$  to  $+40^{\circ}$ . 2-section upholstery 70mm thick, leatherette cover acc. to SCHMITZ colour chart. Powder coated chassis in white aluminium RAL 9006, grey aluminium RAL 9007 or grey white RAL 9002. Silver-grey plastic chassis cover. With 4 levelling feet. Safe working load: 225 kg

### 211.4110.0 Examination and treatment couch<sup>123</sup>

800 mm wide, otherwise like 211.4100.0.



211.4200.0 + 200.0250.0 + 200.0400.0

### 211.4200.0 Examination and treatment couch<sup>123</sup>

electromotive height adjustment from 500 to 920 mm by means of foot control unit, with safety ON/OFF switch to turn control element on or off, 2.40 m connecting cable with mains plug, otherwise like 211.4100.0. Upon request, hand control instead of foot control unit.

### 211.4210.0 Examination and treatment couch<sup>123</sup>

800 mm wide, otherwise like 211.4200.0.

- <sup>1</sup> 4 colours for leatherette upholstery in "classic" quality without surcharge.
- <sup>2</sup> 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.
- <sup>3</sup> Standard powder-coating colours for frame: grey-white 9002, white aluminium RAL 9006 or grey aluminium RAL 9007. Please indicate desired colour upon order placement.

varimed® C€<sup>MD</sup>

## Height adjustable examination and treatment couches

Model No. Description



**200.0250.0 Mobility**

4 closed double castors Ø 100 mm,  
2 of which el. conductive, central lock per axis,  
1 castor with directional lock.  
(Height adjustment range remains the same.)

**Surcharge**

**For mobility for retrofitting,  
please contact our service team.**



**200.0400.0 Side guards, pair**

length: 990 mm, folding down to below top of  
mattress, powder coated white aluminium RAL 9006.



**200.0350.0 Holder for foot control**

for Examination and Treatment couches 211.4200.0  
and 211.4210.0, of steel, powder coated white  
aluminium RAL 9006, for mounting at the long side  
of the couch (retrofittable)

**200.0930.0 Paper roll holder**

(w/o ill.) for paper roll up to 590 mm wide.

**200.1030.0 Paper roll**

(w/o ill.) 500 mm wide

**200.0410.0 Edge protection**

(w/o ill.) of plastic, grey, at all four upholstery corners

**200.0500.0 Head section with nose outlet**

(w/o ill.)

**Surcharge**

**200.0510.0 Filler for nose outlet**

(w/o ill.)

**Surcharge**

Model No.	Description	
<b>200.4402.0</b> (w/o ill.)	<b>Side rail (piece)</b> for Examination and Treatment couches 211.4100.0, 211.4110.0, 211.4200.0 and 211.4210.0, mounting either side, 25 × 10 × 750 mm, max. charge 30 kg, cannot be combined with side guards 200.0400.0. for 211.4100.0, 211.4110.0, 211.4200.0 and 211.4210.0.	
<b>200.0251.0</b> (w/o ill.)	<b>ECG shield</b> with one potential equalization sockets each LHS and RHS as well as screening plate to reduce disturbance, suitable for 650 mm wide couches.	
<b>200.0252.0</b> (w/o ill.)	<b>ECG shield</b> for 800 mm wide couches, otherwise like 200.0251.0.	
<b>200.6530.0</b> (w/o ill.)	<b>Surcharge for coloured chassis</b> acc. to colour chart.	<b>Surcharge</b>



222.2760.0

222.2810.0

222.1550.0

222.1570.0

222.2760.0 upholstery milky coffee 102  
222.2810.0 upholstery raspberry red 103  
222.1550.0 and 222.1570.0 upholstery caribbean blue 82



221.1860.0  
upholstery classic silver grey 55, powder coating grey white RAL 9002




242.2351.0



## varimed® Swivel stools, Examination seats, Medical furniture

<sup>1</sup> classic  
<sup>2</sup> comfort

Model No. Description




### 222.2810.0 Swivel stool<sup>12</sup>

5-arm plastic base, black, chromium-plated column, friction braked double castors for smooth/hard floors, height adjustment by means of "anatomic-lift" gas spring, with hand ring release, chromium-plated column, leatherette upholstered seat. Please specify double castors for use on carpeting upon order, if required. Safe working load: 135 kg.

200.8000.0 Powder coated base according to colour chart  
chromium-plated column.

Surcharge



### 222.2760.0 Swivel stool<sup>12</sup>



5-arm polished aluminium base, chromium-plated column, friction-braked double castors for smooth/hard floors, height adjustment by means of "anatomic-lift" gas spring, with hand ring release, leatherette upholstered seat. Please specify double castors for use on carpeting upon order, if required. Safe working load: 135 kg.

<sup>1</sup> 4 colours for leatherette upholstery in "classic" quality without surcharge.

<sup>2</sup> 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.

## varimed® Swivel stools, Examination seats, Medical furniture

<sup>1</sup> classic  
<sup>2</sup> comfort



Model No. Description

**222.1550.0 Examination stool<sup>12</sup>**  
specially designed upholstered seat with leatherette cover acc. to colour chart, height adjustment from 490 to 690 mm by means of gas spring, with foot ring release. 5-arm base powder coated in grey aluminium RAL 9007, with antistatic plastic castors Ø 50 mm for smooth/hard floors, with perimetral black bumper protection.  
Safe working load: 135 kg.


**222.1570.0 Examination stool<sup>12</sup>**  
specially designed upholstered seat with leatherette cover acc. to colour chart, height adjustment from 490 to 690 mm by means of gas spring, with foot ring release, foam rubber backrest with leatherette cover, 5-arm base powder coated in grey aluminium RAL 9007, with antistatic plastic castors Ø 50 mm for smooth/hard floors, with perimetral black bumper protection.  
Safe working load: 135 kg.

**200.0011.0 Brake loaded castors**  
braked when loaded, released when you come off the stool. For examination seats 222.1550.0 and 222.1570.0. **Surcharge**

<sup>1</sup> 4 colours for leatherette upholstery in "classic" quality without surcharge.

<sup>2</sup> 15 colours for leatherette upholstery in "comfort" quality (surcharge). Please see page 40.

<sup>1</sup> classic  
<sup>2</sup> comfort





Model No. Description

**221.1860.0 Examination/Blood sampling chair<sup>123</sup> €<sup>€MD</sup>**  
base of rectangular tube, powder coated in RAL 9002, RAL 9006 or RAL 9007, seat, backrest and armrests upholstered in leatherette classic black 57. Unassembled delivery. Safe working load: 225 kg.

**200.6610.0 Special coloured finish of frame**  
acc. to colour chart.

**Surcharge**




**242.2351.0 Drip-feed stand €<sup>€MD</sup>**  
 5-arm base of 18/10 stainless steel frame, drip glass, double plastic castors diam. 50 mm, 2 of which lockable, antistatic, upper section with 4 18/10 stainless steel hooks, drip glass, height adjustment from 1,350 to 2,150 mm, one-hand action. footprint 635 mm  
Safe working load: 2 kg (per hook).

- <sup>1</sup> 3 further colours for leatherette upholstery in “classic” quality without surcharge.
- <sup>2</sup> 15 colours for leatherette upholstery in “comfort” quality (surcharge). Please see page 40.
- <sup>3</sup> Standard powder-coating colours for frame: grey-white 9002, white aluminium RAL 9006 or grey aluminium RAL 9007. Please indicate desired colour upon order placement.



medi-net general purpose table 233.6060.0,  
frame grey aluminium RAL 9007

## varimed® medi-net general purpose table



Model No. Description

**233.6060.0 medi-net general purpose table<sup>3</sup>**  
 rectangular steel tube frame, collapsible, powder coated grey white RAL 9002, white aluminium RAL 9006 or grey aluminium RAL 9007, drawer with pre-assembled black galvanized ball bearing precision runners, table top, shelf and drawer front in pearl white 60 plastic veneers, handle for drawer front powder coated in grey white RAL 9002, white aluminium RAL 9006 or grey aluminium RAL 9007, 4 double castors diam. 75 mm.  
 The medi-net table 233.6060.0 can be fitted with 2 more drawers, in addition to the one standard drawer.

Unassembled delivery.

### medi-net

#### Safe working load:

Totally: up to 50 kg

Drawer: up to 10 kg

Top/shelf: up to 20 kg

**201.3000.0 1 additional drawer  
with handle**

**200.6600.0 Coloured frame acc. to colour chart**

**Surcharge**

**200.6650.0 Coloured handle acc. to colour chart**

**Surcharge**

<sup>3</sup> Standard powder-coating colours for frame:  
 9002, RAL 9006 or RAL 9007.  
 Please indicate desired colour upon order placement.



medi-net general purpose table 233.6060.0  
with accessories for injecting workplace:  
201.1170.0, 201.1610.0, 201.1700.0, 201.2600.0



medi-net general purpose table 233.6060.0  
with accessories for ergometry measuring  
workstations 201.1300.0 (w/o equipment)



medi-net general purpose table 233.6060.0 +  
201.0050.0 with accessories for home dialysis  
201.1800.0, 201.1250.0 (w/o equipment)



Cannula and syringe dispenser 201.1610.0



Cannula and syringe dispenser 201.1650.0





Console 201.1800.0



Ampule opener 201.2600.0



“Medi-Müll” disposal set 201.2700.0



1 round SST 18/10 tray 201.1170.0



Infusion pole 201.1250.0





Refuse collector 201.2800.0

**varimed®**

## Accessories for medi-net general purpose table


Model No.	Description	
<b>201.3000.0</b> (w/o ill.)	<b>1 drawer</b> with handle, running slides for subsequent fitting, (internal dimensions: 510 × 413 × 90 mm, W x D x H)	
<b>201.0050.0</b>	<b>2 castors with locking device</b> 75 mm Ø.	<b>Surcharge</b>
<b>201.1300.0</b>	<b>Cable guide</b> with fixing clamp 201.1910.0, 16 mm Ø.	
<b>201.1170.0</b>	<b>1 round SST tray</b> diam. 320 mm, with fixing clamp 201.1900.0, diam. 10 mm	
<b>201.1250.0</b>	<b>Infusion pole <small>CE</small></b> 2-section, with 4 hooks and fixing clamp 201.1920.0, diam. 18 mm, safe working load: 2 kg (per hook)	






Model No. Description


**201.1610.0 Cannula and syringe dispenser**  
grey, first in/first out system, with twin set incl. swab dispenser, requires fixing brackets 201.1700.0




**201.1640.0 Cannula and syringe dispenser**  
grey, first in/first out system, without right twin set, requires fixing brackets 201.1700.0




**201.1650.0 Cannula and syringe dispenser**  
grey, first in/last out system, with twin set incl. swab dispenser, requires fixing brackets 201.1700.0



**201.1660.0 Cannula and syringe dispenser**  
grey, first in/last out system, without twin set, requires fixing brackets 201.1700.0



**201.1700.0 Fixing bracket, 1 pair**  
for use with 201.1610.0, 201.1640.0, 201.1650.0 and 201.1660.0



**varimed®**

## Accessories for medi-net general purpose table

Model No. Description

---

**201.1800.0 Console**  
SST 18/10 with fixing brackets

---

**201.1900.0 Fixing clamp**  
(w/o ill.) 10 mm Ø


---

**201.1910.0 Fixing clamp**  
(w/o ill.) 16 mm Ø


---

**201.1920.0 Fixing clamp**  
(w/o ill.) 18 mm Ø

Model No. Description




**201.2600.0** Ampule opener  
"Cupfix"



**201.2700.0** "Medi-Müll"  
disposal set 1,5 l, 20 boxes made of polypropylene  
with red lid, with holder

**201.2710.0** "Medi-Müll" spare boxes  
(w/o ill.) 1,5 l, carton of 100 boxes, without holder

**201.2720.0** Holder  
(w/o ill.) as 201.2700.0, however without boxes.  
(Holders for different systems upon request)



**201.2800.0** Refuse collector  
with roll of bags, powder coated grey white

**201.2850.0** Roll of spare bags  
(w/o ill.) for 201.2800.0

## **varimed®**

### **Model and material description**

#### **Couches**

Tubular steel frame 40 × 20 mm, pair of feet 30 × 30 mm, powder-coated or chrome-plated. With levelling feet to compensate for floor unevenness. Head section with 2 adjustment mechanisms, 100% metal.

#### **medi-net**

Tubular steel frame 20 × 20 mm, plastic enamel coated.

#### **Castors**

All castors are non-marking and resistant to abrasions.

#### **Drawers**

Drawers run on precision ball-bearing slides.

#### **Handles**

of ergonomic design, distance of holes: 192 mm. The bowed bar handles are chemical resistant, easy to clean and to disinfect.

#### **Powder coated finish**

of high quality, multilayer gloss paint of extreme knock and scratch resistance. Easy to clean. Resistant to climatic conditions and disinfectants.

#### **Standard finish of tubular steel frames**

e.g. for couches: grey white RAL 9002, white aluminium RAL 9006, grey aluminium RAL 9007 or chromium-plated. Please refer to the different model descriptions.

#### **Frames of couches, medi-net trolleys or stool bases**

are also available in other colour acc. to colour chart at extra charge.

#### **Further special powder coatings**

acc. to RAL colour chart on request!



**Upholstery**

Foamed material on a support plate, with first-class, hard-wearing leatherette covers in two qualities.

**Leatherette comfort (available for an additional charge): 15 colour tones**

Wear-resisting material, with high resilience and flexibility providing excellent sitting and lying comfort. Grained surface.

**Leatherette classic (available at no extra charge): 4 colour tones**

Hard-wearing, flexible material. Smooth surface.

**We recommend:****Accessories for upholstery cleaning**

Model No.	Description
<b>2019878</b> (w/o ill.)	<b>Special cleaner for leatherette</b> 500 ml spray bottle (packaging unit 3 bottles)
<b>2026973</b> (w/o ill.)	<b>Cleaning set for leatherette</b> 1 × 500 ml spray bottle and a brush

Please refer to our separate operating manual for upholsteries for detailed information on cleaning and disinfection of upholsteries.

**Leatherette comfort**


gold yellow 86



agave green 109




gold orange 106



milky coffee 102




red carmine 92



sand beige 101




raspberry red 103



stone grey 96




ice blue 108



shale grey 97





caribbean blue 82



titan grey 91



atlantic blue 83

blossom white  
84lavender  
purple 107**Leatherette classic**

white 53



silver grey 55




ocean 52



black 57

**Plastic  
elements**

silver-grey

grey white  
RAL 9002pure white  
RAL 9010alabaster grey  
RAL 0007500<sup>5</sup>light grey  
RAL 7035anthracite grey  
RAL 7016jet black  
RAL 9005white aluminium  
RAL 9006grey aluminium  
RAL 9007colza yellow  
RAL 1009050<sup>5</sup>signal yellow  
RAL 1003pastel orange  
RAL 2003light ivory  
RAL 1015sahara beige  
RAL 0607030<sup>5</sup>antique pink  
RAL 3014heather violet  
RAL 4003vermilion  
RAL 2002light green  
RAL 6027lemon green  
RAL 1008060<sup>5</sup>leaf green  
RAL 1206050<sup>5</sup>white blue  
RAL 2408015<sup>5</sup>light blue  
RAL 2606030<sup>5</sup>sapphire blue  
RAL 5003

antimicrobial  
These colours have antimicrobial properties

## Colours for upholsteries and plastic enamel coatings

### Discover new, colourful worlds

With a total of 19 upholstery colours and 22 colours for powder-coated finish, you have many options to create colourful highlights in your surgery.

Whether muted or cheerful, classic or modern, elegant or extravagant – with our tool, you will certainly find the right colour for you in no time at all.

We would be delighted to advise you – please feel free to contact us.



Be inspired:  
[www.schmitz-soehne.com](http://www.schmitz-soehne.com)



## Contact

SCHMITZ u. Söhne  
GmbH & Co. KG  
Zum Ostefeld 29  
58739 Wickede (Ruhr)

P. O. Box 14 61  
58734 Wickede (Ruhr)  
Germany

T +49 (0)2377 84 0  
F +49 (0)2377 84 162  
zentrale@schmitz-soehne.de

Exports department:  
T +49 (0)2377 84 0  
F +49 (0)2377 84 135  
export@schmitz-soehne.de

Technical service, hotline:  
T +49 (0)2377 84 549  
F +49 (0)2377 84 210  
service@schmitz-soehne.de

Information on  
our establishments  
worldwide can be found  
on our website:

[www.schmitz-soehne.com](http://www.schmitz-soehne.com)



All devices marked with the CE-MD symbol, in combination with our accessories, meet the requirements of Medical Device Regulation 2017/745/EU (MDR) as well as the Medical Device Amendment Act and are risk class 1 medical devices.



All devices marked with the simple CE symbol meet the requirements of low-voltage directive 2014/35/EU as well as the regulation on electrical equipment.

Please also refer to the instructions for use, for more details about the intended use!



\* This product is an article of merchandise, sold under OEM brand name.

Our goods and services are to be understood ex works. Deliveries are effected via our authorized partners or branch offices. The company reserves the right to alter the design, construction, dimensions and finish of its products without notice. All relevant dimensions and weights are to be understood as approximate information/data. Due to technical reasons colour deviations are possible. The articles shown in this catalogue are not authorized for sale in the United States.

Nominal charge 3.00 Euro.