

---

## UniLeg leg supports

Flexible and easy to use, protective for the patients.



- Safe and comfortable leg positioning, foot and lower leg are comfortably padded.
- Extra well-padded peroneal protection shields the sensitive nerves at the knee from compression damage.
- The leg holders can be adjusted conveniently and quickly to any desired position by means of gas spring support.
- Precise attachment to any standard rail by thanks to an integrated clamp standard rail, no matter if examination chair or operating table.
- Made by SCHMITZ - another highlight from our own development, designed and produced in-house.
- Service-friendly through simplified access to the joint.


**101.0369.0**  
**UniLeg leg supports (pair)**

---

**With peroneal nerve protection**, height-adjustable and spreadable from foot end, with gas spring support, plastic covering with pad for foot and lower leg, adjustable in all positions, with integrated side rail mount. Safe working load: 182 kg.

**101.0678.0**  
**UniLeg leg supports (pair)**

Same as model 01.0369.0, but **without peroneal protection**.



SCHMITZ

A SENSE OF CONFIDENCE