

LED Carla (25°)

LiD
Eclairage Lighting 

An examination lamp with a modern and refined design. A compact "water-drop" shaped head enhances ease of use in multiple situations whilst reducing hindrance.

The flexible and silent gooseneck provides for excellent handling, and precise and steady positioning. Very uniform diffuse (25°) lighting.

Ideal as an examination lamp for nurses use, or for any medical speciality and general practice



Carla 65
On/off Switch

● Medical device class 1


Light intensity **420 Lumens***

Colour of temperature **4500°K* (neutral white)**

Colour rendering index **CRI 95/100**

Electric consumption **4.2 Watts**

- LED life span: 50 000 h minimum
- Electrical protection: insulation class II
- Standards: IEC 60 601 -1/ -2-41
- Power cable length: 3 meters
- Suggested colors: white and black
- Head dimensions: 6 x 10 cm
- Fittings: 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13

Lux*	Beam 25°	Distance	Illuminated diameter
20000		30 cm	15 cm
7000		50 cm	20 cm
1700		100 cm	30 cm

* Values accurate within +/- 10 %

CODES		
	On/off Switch	Contactless switch
Flexible 65 cm	LEDB0650	LEDB0650-ISC
Flexible 100 cm	LEDB0900	LEDB0900-ISC
Flexible 123 cm	LEDB1200	LEDB1200-ISC

